

SURGEON - PATRICK WEINRAUCH**PROCEDURE - TOTAL HIP ARTHROPLASTY - ANTERIOR APPROACH**

POSITION - SUPINE

<p style="text-align: center;"><u>GLOVES</u></p> PATRICK - 2 x 8 BLUE & PURPLE SHARON - 3 x 7 (1 X BLUE & 2 X PURPLE) CARMEN - 2 X 7 BLUE & PURPLE	<p style="text-align: center;"><u>DRAPES</u></p> UNIVERSAL PACK, U-BAR PACK, BILATERAL EXTREMITY DRAPE, LARGE DRAPE, 1/2 DRAPE, MAYO COVER, LEGGINGS, LARGE STOCKINETTE X 4, LARGE OPSITE
<p style="text-align: center;"><u>GOWNS</u></p> 1/2 DRAPE, XL MICROCOOL X 2, L MICROCOOL X 1 STRYKER HOODS X 3	FOR BILATERAL PROCEDURE SEE OVER PAGE . .
<p style="text-align: center;"><u>BOWLS</u></p> MAJOR BOWL SET, JUG	<p style="text-align: center;"><u>SKIN PREP</u></p> ALCOHOLIC CHLORHEXADINE
<p style="text-align: center;"><u>EXTRAS</u></p> PREP SPONGE x 2 PURPLE ALEXIS RAYTEC X 1 COLLOTAMP 10 X 10 MEDIUM SPONGES X 2 CEPHAZOLIN 1g (WASH) 10cm GAUZE RIBBONS X 2 SUCTION TUBING NEEDLE MAT LIGHT HANDLE 14g ORANGE JELCO 3ml SYRINGE 50ml SYRINGE X 2 DIATHERMY & MEDIUM INSULATED TIP KAM-VAC SUCTION BULB SYRINGE DRAPE TAPE	<p style="text-align: center;"><u>INSTRUMENTS</u></p> MAJOR ORTHO TRAY TRS DRILL AND BATTERY ANTERIOR HIP EXTRAS LARGE MALLET NORFOLK NORRICH MAGNETIC MAT LONG, WIDE LANGENBECKS HIP CAPSULE SPREADERS (CAN USE LINK SHOULDER TRAY) OT 8 HARD LIGHT HANDLE LOAN SETS
<p style="text-align: center;"><u>SUTURES AND INFILTRATION</u></p> 0 MAXON LOOP (GMM563L) x 2 0 V-LOC (N0326) X 1 3/0 VLOC 90 (M0024) X 1 LIQUIBAND NAROPIN 0.2% (200mls), ADRENALINE, TORADOL +/- MORPHINE (CHECK EXACT COCKTAIL W ANAESTHETIST)	<p style="text-align: center;"><u>DRESSINGS</u></p> MEPILEX BORDER
<p style="text-align: center;"><u>DISPOSABLES</u></p> SAW BLADE (REF- 05.002.003S)	<p style="text-align: center;"><u>EQUIPMENT</u></p> STERIS TABLE, FULL LENGTH JELLY MAT NEPTUNE SUCTION, DIATHERMY (CUT 70, COAG 70) YELLOW FINS ORT LATERAL HIP SUPPORTS STERIS KNUCKLES X 4 FLUID TROLLEY BLACK ANTI FATIGUE MATS

SURGEON - PATRICK WEINRAUCH

PROCEDURE - TOTAL HIP ARTHROPLASTY - ANTERIOR APPROACH

POSITION - SUPINE

DRAPING ORDER

- | | |
|--|---|
| 1. LEGGING ON YELLOW FIN | 11. OPSITE |
| 2. BOTTOM DRAPE | 12. GAUZE RIBBONS - 2 X 1/2 LENGTHS ON OPERATIVE
SIDE SUPPORTS, 1 X FULL LENGTH FOLDED IN HALF AND
PLACED ON OPPOSITE POST, 1 X ROBERTS AND 2 X
HARRISONS PLACED ON EACH RIBBON AND SECURED TO
DRAPES |
| 3. 2 X STOCKINETTE ON NON OPERATIVE LEG | 13. DIATHERMY/SUCTION PLACED DOWN LENGTH OF PT. |
| 4. SPLIT DRAPE | 14. TOP DRAPE ACROSS PT TO SECURE LINES |
| 5. TOP DRAPE FROM U-BAR PACK | 15. 2 X SIDE DRAPE SECURING TOP DRAPE AND RUNNING
THE LENGTH OF THE PATIENT |
| 6. 2 X STOCKINETTE ON OPERATIVE LEG | |
| 7. OPSITE TO OPERATIVE AREA | |
| 8. BILATERAL EXTREMITY DRAPE | |
| 9. CUT BILATERAL DRAPE OVER OPERATIVE AREA | |
| 10. SECOND LEGGING | |

FOR BILATERAL ADD:

- | | |
|------------------------|----------------------|
| U-DRAPE X 1 | COLLATAMP 10 X 10 |
| 22 BLADE X 1 | SUTURES AND DESSINGS |
| LARGE OPSITE X 1 | PURPLE ALEXIS |
| 10cm GAUZE RIBBONS X 2 | |

GENERAL TIPS

- * INTRA-OP LOCAL INFILTRATION COCKTAIL - 150mls NAROPIN 0.2% +/- ADRENALINE +/- TORADOL +/- MORPHINE
CHECK WITH ANAETHETIST BEFORE EACH CASE
- * 1G CEPHAZOLIN IN WASH
- * FOR BILATERAL PROCEDURES SPLIT LOCAL COCKTAIL BETWEEN BOTH SIDES - 75mls EACH SIDE
- * DIATHERMY AND SUCTION ON OPPOSITE SIDE FROM ANAESTHETICS REGARDLESS OF SIDE BEING OPERATED ON