

SURGEON - PATRICK WEINRAUCH PROCEDURE - REVISION TOTAL HIP ARTHROPLASTY POSITION - LATERAL	
<p align="center"><u>GLOVES</u></p> PATRICK - 2 x 8 BLUE & PURPLE SHARON - 3 x 7 (1 X BLUE, 2 X PURPLE) CARMEN - 2 X 7 (BLUE & PURPLE)	<p align="center"><u>DRAPES</u></p> HIP PACK, LARGE DRAPE x 2, HALF DRAPE, OPSITE x 2 IMPERVIOUS SPLIT DRAPES x 1, LARGE STOCKINETTE MAYO COVER x 2, BACK TABLE COVER x 2
<p align="center"><u>GOWNS</u></p> 1/2 DRAPE, XL MICROCOOL X 2, L MICROCOOL X 1 STRYKER HOODS X 2	
<p align="center"><u>BOWLS</u></p> MAJOR BOWL SET, COLLAPSIBLE BOWL	<p align="center"><u>SKIN PREP</u></p> ALCOHOLIC CHLORHEXADINE
<p align="center"><u>EXTRAS</u></p> PREP SPONGE x 2 MEDIUM SPONGE x 2 LARGE SPONGES x 2 SUCTION TUBING x 2 NEEDLE MAT LIGHT HANDLE x 2 22 BLADE x 2 50ml SYRINGE x 2 3 ml SYRINGE 18G BLUNT NEEDLE 14G ORANGE JELCO DIATHERMY & MEDIUM INSULATED TIP SURGILAV KAM-VAC SUCTION STERILE SPECIMEN POTS X 3 PATHOLOGY SWABS X 3 5 pk BLUE TOWELS <u>HAVE AVAILABLE:</u> CEMENT BOWLS (DNO) ANTIBIOTIC SIMPLE CEMENT (DNO) LARGE IN/OUT CATHETER (DNO) 5cm & 10cm RIBBON GAUZE PACK (DNO) 10fg BELLOVAC DRAIN X 2	<p align="center"><u>INSTRUMENTS</u></p> MAJOR ORTHO TRAY TRS DRILL AND BATTERY HEAVY MALLET POSTERIOR HIP EXTRAS TRAY NORFOLK NORRICH x 2 CHARNLEY BOW RETRACTOR LOAN INSTRUMENTS <u>HAVE AVAILABLE:</u> DALL MILES (DNO) SURGIATOME AND CABLE (DNO) CEMENT COMPRESSORS (DNO) SWISS BONE PUNCHES (DNO) BONE MILL (DNO) LARGE FRAGMENT SCREWS (DNO) HIBBS RETRACTORS (DNO) LARGE BONE NIBBLERS (DNO) WIRE CUTTERS (DNO) WIRE PLIERS/MUTIGRIPS (DNO) HAYGROVE BONE HOLDERS (DNO) FRABRUGE BONE HOLDERS (DNO) LARGE BONE REDUCTION FORCEPS (DNO) NARROW HOHMANS (DNO) BROAD HOHMANS (DNO)

SURGEON - PATRICK WEINRAUCH

PROCEDURE - REVISION TOTAL HIP ARTHROPLASTY

POSITION - LATERAL

SUTURES AND INFILTRATION

0 MAXON LOOP (6325-63) x 3
0 VICRYL (346) x 4
3/0 VLOC 180 (19mm NEEDLE)
SKIN STAPLES x 2
DERMABOND

NAROPIN 0.2% (200mls), ADRENALINE +/- TORADOL +/-
MORPINE - CHECK WITH ANAESTHETIST

DRESSINGS

LARGE PRIMAPORE x 4
COMBINE 20 x 30cm x 2
HYPERFIX x 2

DISPOSABLES

SAW BLADE (REF- 519.113S)
ANTIBIOTIC SIMPLEX (DNO)
BEELINE PAIN PUMP (250mls)

EQUIPMENT

STERIS TABLE
CELL SALVAGE
NEPTUNE SUCTION
DIATHERMY (CUT 70, COAG 100)
ORT LATERAL SUPPORTS
STERIS KNUCKLES
FLUID TROLLEY
BLACK ANTIFATIGUE MATS

GENERAL TIPS

- * 1 AMP CEPHAZOLIN TO EACH 2L BAG OF NORMAL SALINE
- * CELL SALVAGE WILL BE USED FOR ALL REVISION PROCEDURES
- * NEED THE TROLLEY WITH BRAKES ON THE FEET FOR THE BONE MILL, REINFORCE WITH SANDBAGS
- * DIATHERMY AND SUCTION ON PATRICK'S SIDE (PATIENT'S BACK)